





# Acrylic anhydride

# Methacrylic anhydride

## Features

- ① High quality at 99% purity or more.
- ② 99.9% purity or more is also available.
- ③ Residual chlorine less than 0.005%
- ④ Free acid less than 0.1%
- ⑤ High stability (Refrigerated storage: 6 months)
- ⑥ Polymerization inhibitors can be customized.


## Physical properties

Acrylic anhydride	Product name	Methacrylic anhydride																																																																																																						
$C_6H_6O_3$	Molecular formula	$C_8H_{10}O_3$																																																																																																						
126.11	Molecular weight	154.17																																																																																																						
	Structural formula																																																																																																							
2051-76-5	CAS No.	760-93-0																																																																																																						
Colorless transparent liquid	External appearance	Colorless transparent liquid																																																																																																						
Irritating odor	Odor	Irritating odor																																																																																																						
84°C/2kPa	Boiling point	90°C/2kPa																																																																																																						
-20°C	Melting point (lit)	<-25°C																																																																																																						
1.094/20°C	Specific gravity (lit)	1.032/20°C																																																																																																						
99% or higher	Purity (GC)	99% or higher																																																																																																						
MEHQ 100ppm	Additives	MEHQ 100ppm																																																																																																						
Representative chromatography																																																																																																								
 <table border="1"> <thead> <tr> <th>ピーク番号</th> <th>保持時間</th> <th>高さ</th> <th>面積</th> <th>面積%</th> <th>化合物名</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>5.806</td> <td>1714</td> <td>13848</td> <td>0.0120</td> <td></td> </tr> <tr> <td>2</td> <td>8.025</td> <td>3714</td> <td>30240</td> <td>0.0261</td> <td></td> </tr> <tr> <td>3</td> <td>10.897</td> <td>4285029</td> <td>115521390</td> <td>99.8978</td> <td>アクリル酸無水物</td> </tr> <tr> <td>4</td> <td>12.050</td> <td>214</td> <td>1582</td> <td>0.0014</td> <td></td> </tr> <tr> <td>5</td> <td>15.454</td> <td>19482</td> <td>63149</td> <td>0.0546</td> <td></td> </tr> <tr> <td>6</td> <td>19.234</td> <td>395</td> <td>1926</td> <td>0.0017</td> <td></td> </tr> <tr> <td>7</td> <td>19.551</td> <td>694</td> <td>7475</td> <td>0.0065</td> <td></td> </tr> <tr> <td>合計</td> <td></td> <td>4311242</td> <td>115639610</td> <td>100.0000</td> <td></td> </tr> </tbody> </table>	ピーク番号	保持時間	高さ	面積	面積%	化合物名	1	5.806	1714	13848	0.0120		2	8.025	3714	30240	0.0261		3	10.897	4285029	115521390	99.8978	アクリル酸無水物	4	12.050	214	1582	0.0014		5	15.454	19482	63149	0.0546		6	19.234	395	1926	0.0017		7	19.551	694	7475	0.0065		合計		4311242	115639610	100.0000			 <table border="1"> <thead> <tr> <th>ピーク番号</th> <th>保持時間</th> <th>高さ</th> <th>面積</th> <th>面積%</th> <th>化合物名</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>5.502</td> <td>398</td> <td>1228</td> <td>0.0010</td> <td></td> </tr> <tr> <td>2</td> <td>7.188</td> <td>5749</td> <td>49570</td> <td>0.0398</td> <td></td> </tr> <tr> <td>3</td> <td>11.529</td> <td>1249</td> <td>14667</td> <td>0.0118</td> <td></td> </tr> <tr> <td>4</td> <td>12.506</td> <td>1063</td> <td>13104</td> <td>0.0105</td> <td></td> </tr> <tr> <td>5</td> <td>13.660</td> <td>5943916</td> <td>124457432</td> <td>99.9341</td> <td>メタクリル酸無水物</td> </tr> <tr> <td>6</td> <td>20.124</td> <td>790</td> <td>3508</td> <td>0.0028</td> <td></td> </tr> <tr> <td>合計</td> <td></td> <td>5953165</td> <td>124539509</td> <td>100.0000</td> <td></td> </tr> </tbody> </table>	ピーク番号	保持時間	高さ	面積	面積%	化合物名	1	5.502	398	1228	0.0010		2	7.188	5749	49570	0.0398		3	11.529	1249	14667	0.0118		4	12.506	1063	13104	0.0105		5	13.660	5943916	124457432	99.9341	メタクリル酸無水物	6	20.124	790	3508	0.0028		合計		5953165	124539509	100.0000	
ピーク番号	保持時間	高さ	面積	面積%	化合物名																																																																																																			
1	5.806	1714	13848	0.0120																																																																																																				
2	8.025	3714	30240	0.0261																																																																																																				
3	10.897	4285029	115521390	99.8978	アクリル酸無水物																																																																																																			
4	12.050	214	1582	0.0014																																																																																																				
5	15.454	19482	63149	0.0546																																																																																																				
6	19.234	395	1926	0.0017																																																																																																				
7	19.551	694	7475	0.0065																																																																																																				
合計		4311242	115639610	100.0000																																																																																																				
ピーク番号	保持時間	高さ	面積	面積%	化合物名																																																																																																			
1	5.502	398	1228	0.0010																																																																																																				
2	7.188	5749	49570	0.0398																																																																																																				
3	11.529	1249	14667	0.0118																																																																																																				
4	12.506	1063	13104	0.0105																																																																																																				
5	13.660	5943916	124457432	99.9341	メタクリル酸無水物																																																																																																			
6	20.124	790	3508	0.0028																																																																																																				
合計		5953165	124539509	100.0000																																																																																																				



## As various reaction reagents

- (Meth) acrylic anhydride reacts with functional group compounds to form derivatives and is used for various applications.



### Liquid crystal materials (liquid crystal sealants, UV curable resins)

- Polyfunctional esters synthesized by reaction with epoxy resin are used as liquid crystal sealants, UV curable resins, etc.

### Photoresist

- Alicyclic esters synthesized by reaction with alcohols, etc. are used for photoresist applications.

### Paints, resins, adhesives, etc.

- Amides and esters obtained by reaction with an amino group or hydroxyl group are used in paints, resins, adhesives, etc.

### Optical Lenses

- Esters obtained by reaction with fluorene derivatives and aromatic compounds such as naphthol, are used as a high refractive lens monomers.

### Photoinitiator

- Compounds obtained by reaction with photosensitive phenol derivatives are used as photoinitiators for UV paints and coatings.

### Medical and agrochemicals

- Amide derivatives obtained by reaction with benzylamine are used as intermediates for antiepileptic drugs.
- Esters obtained by reaction with lactone compounds, which have a hydroxyl group, are used in the medical and agrochemical fields.

### Dental materials

- Esters obtained by reaction with a benzoate alcohol derivatives are used as a dental adhesive monomers.

### Fuel cells

- Multifunctional esters synthesized by reaction with epoxy resin are used as separators for fuel cells.

### Water-based sizing agent

- Amide derivatives obtained by reaction with acrylamide-functionalized alkoxy silane are used as water-based sizing agents

### Polymeric gas hydrate formation inhibitor

- Esters obtained by reaction with N-alkylamine are used as building blocks for polymeric gas hydrate inhibitors

### Concrete superplasticizer

- Esters obtained by reaction with polyalkylene glycol compounds with a hydroxyl group are used as concrete superplasticizers.